

# Wheelchair Recliner™

## Description

The Recliner is used in clinics where hand wheelchairs predominate. The wheelchair is moved into position and tilted by the device safely and under motor control. The 'Recliner' can be easily moved to a storage position taking minimum space. It can be easily retrofitted to an existing surgery sharing the normal chair's facilities. The 'Wheelchair Recliner™' can be moved between surgeries.

The unit has been introduced to respond to problems experienced in existing surgeries where refurbishment is not an option, space is limited and financial concerns apply constraints. The 'Recliner' eliminates most problems providing a good level of service.

## Specifications

- For use with any standard wheelchair or hospital chair
- Motorised for stress free operation
- Head positioned correctly for comfortable surgery
- Fully adjustable headrest
- Open design for easy access to patient
- Treatment in seated position
- Lift and wheel away
- Ready for use in a minute
- Control by fixed switches
- Optional radio handset control
- Small storage area 200mm x 1500mm
- Minimum storage area allows full surgery access when not in use
- Stands alongside a conventional dental chair to share the services
- Can be moved between surgeries
- Battery operated with overnight charging
- Ideal for mobile surgery vehicles
- Easily fitted in an existing surgery with no additional services required
- 24v operation
- Safeguarded to the machinery directive
- Complies with the Low voltage directory
- Maximum patient weight 127 kg (20 stone)

## Installation

- Floor mounting plates need to be fitted for each operating position
- Site plans available on request
- On site training

## Finance

Through the First National Bank, Design Specific offers Lease Purchase and Rental arrangements. This significantly changes the funding for a project for the repayments can be taken from a revenue rather than capital budget.

